

## Publication 2005

1. "Pyrene-labeled G-quadruplex Oligonucleotide as a Fluorescence Probe for Potassium Ions Detection in Biological Applications" Satoru Nagatoishi, Takahiko Nojima, Bernard Juskowiak & Shigeori Takenaka, *Angew. Chem. Int. Ed.*, **44**, 5067-5070 (2005).
2. "Synthesis of ferrocenylcarbodiimide as a convenient electro-chemically active labeling reagent for nucleic acids" K. Mukumoto, T. Nojima, & S. Takenaka, *Tetrahedron*, **61**, 11705-11715 (2005).
3. "Electrochemical Telomerase Assay with Ferrocenyl naphthalene Diimide as a Tetraplex DNA-Specific Binder" S. Sato, H. Kondo, T. Nojima, & S. Takenaka, *Analytical Chemistry*, **77**, 7304-7309 (2005).
4. "Supramolecular complex formation by beta-cyclodextrin and ferrocenyl naphthalene diimide-intercalated double stranded DNA and improved electrochemical gene detection" Shinobu Sato, Takahiko Nojima, Michinori Waki & Shigeori Takenaka, *Molecules*, **10**, 693-707 (2005).
5. "Supramolecular Assembly of Fullerene Derivatives in the Absence or Presence of Double Stranded DNA in Water" Takahiko Nojima, Kenichi Yamashita, Taizo Hatta, Otohiko Tsuge, Toshio Iwase, Naoko Makita, Kenichi Yoshikawa, Satoshi Fujii & Shigeori Takenaka, *Bunksei Kagaku*, **54**(6), 449-454 (2005).
6. "Genotyping of human lipoprotein lipase gene by using ferrocenyl naphthalene diimide-based electrochemical hybridization assay (FND-EHA)" T. Nojima, K. Yamashita, A. Takagi, M. Takagi, Y. Ikeda, H. Kondo & S. Takenaka, *Analytical Sciences*, **21**, 1437-1441 (2005).
7. "Immobilization of monoclonal antibody for the carcinoembryonic antigen (CEA) on a gold surface and its interaction analysis with CEA using Fourier transform infrared reflection-absorption spectroscopy (FT-IR RAS)" Keiichi Ohtsuka, Masahide Kuroki, Takahiko Nojima, Michinori Waki & Shigeori Takenaka, *Analytical Sciences*, **12**, 215-218 (2005).

8. "Investigation of ferrocenyl carbodiimide (FCDI) in the modification reaction of nucleic acids" Kosuke Mukumoto, Takahiko Nojima, & Shigeori Takenaka, *Nucleic Acids Symposium Series*, **49**, 231-232 (2005).
9. "Electrochemical detection of telomeric quadruplex DNA using ferrocenyl naphthalene diimide" Shinobu Sato, Hiroki Kondo, Takahiko Nojima, & Shigeori Takenaka, *Nucleic Acids Symposium Series*, **49**, 237-238 (2005).
10. "Electrochemical Detection of a Disease-specific Isoform of the Prion Protein (PrP<sup>sc</sup>)" Keiichi Ohtsuka, Katsumi Doh-ura, Michinori Waki, Takahiko Nojima & Shigeori Takenaka, *Peptide Science 2004*, 305-306 (2005).
11. "Fluorescence enhancement of bis-acridine orange peptide, BAO, upon binding to double stranded DNA" Keiji Mizuki, Yutaka Sakakibara, Hiroyuki Ueyama, Takahiko Nojima, Michinori Waki & Shigeori Takenaka, *Organic & Biomolecular Chemistry*, **3**, 578-580 (2005).
12. "Total Analysis Systems, Micro" Shigeori Takenaka, in *Encyclopedia of Molecular Biology and Molecular Medicine* (R. A. MEYERS, Ed.), WILEY-VCH, Weinheim, pp. 421-436 (2005).
13. "Threading intercalators as redox indicators" Shigeori Takenaka, in *Electrochemistry of nucleic acids and proteins. Towards electrochemical sensors for genomics and proteomics* (E. Palecek, F. Scheller, J. Wang, Eds.), Elsevier (2005).